



# Manual Wheelchairs

Manual wheelchairs are available in a number of different sizes, prices and configurations. You will need to consider the following questions to help determine the most suitable model.

## DO YOU NEED TO TRANSPORT THE WHEELCHAIR?

If so, the weight of the wheelchair will be important to the person lifting it in and out of the vehicle.

## STEEL OR ALUMINIUM?

Aluminium wheelchairs are approximately half the weight of the same size steel model. You are also able to remove the wheels and footplates of the alloy models to further reduce the weight.

## WHICH TYPE AND SIZE?

Self propel (24" wheel) are those that can be propelled by the user. Transit versions (12" wheel) must be pushed at all times. Standard sizes available range from 16" wide by 16" deep up to 20" wide by 20" deep. Larger, Bariatric, models are also available upon request. For full details we advise you to talk to one of our trained consultants. They can help you find the right manual wheelchair, at the right price, for your specific needs.

# Manual Wheelchairs

**TRANSIT OR ATTENDANT PROPELLED:** - for those unable to propel themselves or for short transportation of patients over even surfaces. All wheelchairs in this range fold for transportation.

## PRIDE AIRGO NAVIGATOR

Innovative, 2-in-1, walker and transit manual wheelchair in one convenient package. Converts from one to the other in one quick step. The Navigator easily folds up into a small, portable, configuration with carry strap - total weight only 9.5kg. Other features include; tool-free assembly, locking brakes, removable legrests, nylon carry basket, height adjustable handle bars and 113kg weight capacity. Suitable for temporary use only on flat surfaces. Available in various sizes: from 16" to 18" seat width (41-46cm)

Product Code: PMW-NAVIGATOR DVA Code: AP1866116



## PRIDE BRHT TRANSIT

Just released, lightweight transit wheelchair with swing-away, removable legrests, removable armrests, folding back and 12" rear solid tyres. Suitable for temporary use only on flat surfaces for users up to 113kg. Available only in 18" x 16" (46 x 41cm)

Product Code: PMWBRHT



## M4U TRANSIT

Lightweight, economical, aluminium wheelchair. Features include; adjustable and removable swing-away legrests, removable padded armrests, folding back, lockable attendant handbrakes and seat belt. Available in various sizes: from 17" to 22" seat width (43-56cm)

Product Code: MFY-MW19





## Manual Wheelchairs - Transit Continued..

### INVACARE ATLAS LITE TRANSIT

Made from strong, lightweight steel, the Atlas Lite Transit offers full length removable padded armrests, adjustable and removable swing-away legrests, folding back, 12" solid rear tyres and a weight capacity of 120kg. Available in various sizes: from 16" to 19" seat width (41-48cm)

Product Code: *MT-AtlasT+Size?*



### SOMA ECONOMY TRANSIT

Quality, lightweight aluminium alloy manual wheelchair available in various sizes and offering numerous adjustments to suit any individual. Features include; swing-away and removable legrests, flip-back armrest for easy access to desks as well as easy side transfers, adjustable seat height, adjustable front polyurethane castors and pneumatic rear wheels. Available in various sizes: from 16" to 18" seat width (41-46cm)

Product Code: *FL-SomaT+Size ?*



### INVACARE ACTION 2 TRANSIT

Extremely lightweight (just 14.2kg), aluminium wheelchair. Very easy to fold and transport. Features include; swing-in/swing-out legrests with 'quick release' mechanism for maximum manoeuvrability, full or desk length flip-back armrests, 120kg weight capacity, fabric upholstery, 20cm (8") solid front castors and 12.5" pneumatic rear wheels. Available in various sizes: from 16" to 18" seat width (41-46cm)

Product Code: *MT-Action2T+size?*



### KARMA EAGLE TRANSIT II

Top of the range lightweight, strong and durable transit wheelchair offering exceptional comfort. This quality wheelchair has an extensive selection of options and accessories available. Features include; double cross bracing for extra strength and stability, swing back, height adjustable armrests, height adjustable swing-away removable legrests and folding backrest. With a weight capacity of 130kg, 12.5" pneumatic rear wheels and a total weight of between 13.5-14.3kg it represents a highly flexible and robust wheelchair to suit any individual. Available in various sizes: from 16" to 22" seat width (41-56cm)

Product Code: *FL-KarmaT+Size ?*





## Self Propelled

**SELF PROPELLED:** for those wanting to propel the chair themselves or for attendant control especially over uneven ground. All wheelchairs in this range fold for transportation.

### PRIDE XRHT

Economical, steel framed, manual wheelchair primarily suitable for indoor or institutional use due to the heavier framework and design of the tyres. Features include; removable padded armrests, adjustable and removable swing-away legrests, 24" black solid rear tyres, total weight 17kg and a weight capacity of 113kg. Available only in 18" x 16" (46 x 41cm)

*Product Code: PMWXRHT*



### PRIDE QRHT

Offering all the features of the XRHT but with the additional benefit of quick-release rear spoked wheels. The QRHT also offers; removable padded armrests, adjustable and removable swing-away legrests, total weight 18kg and a weight capacity of 113kg. Available only in 18" x 16" (46 x 41cm)

*Product Code: PMWQRHT*



### INVACARE VERANDA

Standard, self-propelled, manual wheelchair made from powder-coated steel. Features include; stylish silver vein frame with matching front riggings, swing-away and removable legrests, padded nylon upholstery, removable padded armrests, total weight of 18kg (w/o legrests) and weight capacity of 113kg. The Veranda is available with the choice of solid or pneumatic 24" tyres and swing-away or padded elevating legrests. Available only in 18" x 16" (46 x 41cm)

*Product Code: MT-Veranda+Size ?*



### INVACARE ATLAS LITE

Compact, lightweight steel, wheelchair with dual-position axles to offer two seat-to-floor heights. Standard features include; swing-away and removable legrests, removable backrest, removable padded armrests, total weight of 16kg (9kg heaviest part) and weight capacity of 120kg. The quick-release rear wheels ensure the Atlas Lite is easy to handle, transport and store. Available in various sizes: from 16" to 19" seat width (41-48cm)

*Product Code: MT-AtlasS+Size ?*





## Self Propelled Continued..

### SOMA ECONOMY 800

Quality, lightweight aluminium alloy manual wheelchair available in various sizes and offering numerous adjustments to suit any individual. Features include; adjustable, swing-away and removable legrests, calf-support strap, flip-back armrest for easy access to desks as well as easy side transfers, adjustable seat height, adjustable front polyurethane castors and pneumatic quick-release rear wheels. With total weights from 14.6kg and a weight capacity of 115kg it represents great value for money. Available in various sizes: from 16" to 18" seat width (41-46cm)

Product Code: FL-SomaS+Size ?



### PRIDE PMW L400

Lightweight, aluminium wheelchair, which folds for easy storage and transportation. It features flip-up, desk-length, armrests, removable swing-away legrests and a height adjustable back. The 24" (61cm) rear wheels are solid and fitted with quick release rear wheels which are set in multi-positional axle plates. It also offers 20cm (8") front solid castors, adjustable rear anti-tip wheels, dracon fabric upholstery and now comes with a deluxe Synergy Solution gel seat cushion for unrivalled comfort and postural support. Available only in 18" x 16" (46 x 41cm)

Product Code: PMWL400



### INVACARE ACTION 2

Extremely lightweight (just 14.2kg), aluminium wheelchair. Very easy to propel or to fold and transport. Features include; swing-in/swing-out legrests with 'quick release' mechanism for maximum manoeuvrability, full or desk length flip-back removable armrests, 120kg weight capacity, fabric upholstery, 20cm (8") solid front castors and 24" (61cm) pneumatic, quick-release rear wheels. Available in various sizes: from 16" to 19" seat width (41-48cm)

Product Code: MT-Action2S+size?



### KARMA S-ERGO RANGE 105,115, 305

The revolutionary-designed S-Ergo 115 offers the unique S-shaped seat that actually contours to the shape of the user to provide pressure relief and prevent slippage. The seat upholstery and removable pads are treated with Aegis Microbe shield to reduce stains and odour, as well as preventing deterioration from bacteria or fungus. The lightweight frame is made from aircraft aluminium and breaks down to just 8kg (total weight 13.1kg). Other features include; streamline oval tubes for added strength, folding back, wide contoured armrests, lockable attendant handbrakes, solid tyres, swing-in/swing-out removable legrests and calf support strap. With a weight capacity of 115kg this stylish wheelchair represents a great all-round package. Available in various sizes: from 16" to 18" seat width (41-46cm) The S-Ergo 105 and 305 models are also available

Product Code: FL-KarmaSES+size? DVA Code: AP17 KM-S+Size?





## Self Propelled Continued..

### KARMA EAGLE

Top of the range lightweight, strong and durable manual wheelchair offering exceptional comfort. The Eagle offers an impressive range of standard features with an extensive selection of personal customising options. Features include; aircraft grade aluminium frame, double cross bracing for extra strength and stability, swing back, height adjustable armrests, Dacron upholstery, height adjustable swing-away removable legrests, weight capacity of 130kg, 24" pneumatic, quick-release rear wheels and a total weight of between 14.4-14.6kg. Available in various sizes: from 16" to 22" seat width (41-56cm)

**Product Code:** FL-KarmaES+Size ? **DVA Code:** AP17 69+Size?



**TILT AND RECLINE CHAIRS:** for those likely to spend extensive periods of time in the chair and requiring high levels of comfort, support and flexibility. All models can be modified to incorporate pressure care cushioning systems.

### RAPTOR TILT & RECLINE

A rigid, steel framed wheelchair with tilt-in-space and reclining back functions. Available only in a transit rear wheel drive model. Features include; in-built adjustability in width and height, adjustable headrest, swing-away legrests (or optional elevating legrests), adjustable armrests and padded, contoured foam upholstery which is removable and washable.

The backrest is height adjustable and may be manually reclined via a lever mechanism and locked in any position from 91 to 108 degrees. The wheelchair features a gas strut tilt-in-space function from 1 to 17 degrees and is operated by an attendant lever.

The seat comes standard with a padded upholstered removable pommel, and padded upholstered removable lateral supports which are height and width adjustable. Seating can be fully customised to suit any user. Weight capacity 160kg and a total weight 42kg, Available in various sizes: from 16" to 17.5" seat width (40-45cm)

**Product Code:** MFY-NW34



### BREEZY RELAX TILT & RECLINE

The Relax is a multi-function tilt and recline (including anterior tilt) wheelchair with a wide range of features and adjustments. It has almost unlimited possibilities to meet the most varied requirements.

It has been designed with functionality in mind with anatomically shaped seat and backrest, adjustable seat height and depth, adjustable height removable armrests with sliding arm pads, elevating detachable legrests with angle adjustable footplates, reclining height adjustable backrest with removable adjustable three part headrest, adjustable wheelbase and frame depth, and quick release transit or self propel rear wheels.

The Dartex upholstery and a combination of high density foam as well as layer of memory foam in the seat ensure a high level of comfort. Seating can be fully customised to suit any user. Weight capacity 140kg and a total weight 38kg. Available in various sizes: from 16" to 20" seat width (41-51cm)

**Product Code:** LD-Relax 309-36

