

PRIDE[®]
MOBILITY



<20 km

Range per Charge



181 kg

Weight Capacity



92.5 cm

Turning Radius



Jazzy[®]
carbon HD

Heavy Duty. Lightweight.

Weighing just 14.5 kg and boasting a 181 kg weight capacity, the Jazzy Carbon HD is a game-changing heavy-duty travel power wheelchair. It folds and locks into position in mere seconds and will keep you comfortable with a padded seat back and extra-wide foot platform.



Visit pridemobility.com.au for more information

(JAZZY-CARBON-HD)



Full front & rear suspension



Quick disconnect joystick with LCD display



Puddle light

RANGE PER CHARGE (MAX) ¹	Up to 20 km
MAXIMUM WEIGHT CAPACITY ³	181 kg
TURNING RADIUS	92.5 cm
MAXIMUM SPEED ¹	6 km/h
LENGTH & WIDTH ²	L: 93.5 cm W: 67 cm
WEIGHT OF HEAVIEST PIECE	14.2 kg
GROUND CLEARANCE	2.5 cm (at motor)
SEAT SIZE ⁴ WIDTH & DEPTH	W: 51 cm x D: 51 cm (20" x 20")
HIGH BACK SEAT-TO-FLOOR RANGE	50 cm fixed

1. Varies with user weight, terrain type, battery amp hour (Ah), battery charge, battery condition and tyre condition.
2. Width of base (not including seat).
3. Maximum weight capacity is reduced to 181 kg on 18" wide seat configuration.

STANDARD ELECTRONICS	YT Control
DRIVE TRAIN	Two motor, rear-wheel drive
BRAKING SYSTEM	Regenerative and electro-mechanical
SUSPENSION TYPE	front and rear
FRONT CASTORS	6.5" solid castors
DRIVE WHEELS	10.6" solid castors
BATTERY REQUIREMENTS	(1) 20Ah, 24V (lithium-ion, 436Wh)
BATTERY CHARGER	2A, off-board
BATTERY WEIGHT	3.4 kg each
COLOUR OPTIONS	■ Black Cherry Carbon

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.