






## THERA-Trainer Tigo – Lower accessories

 A close-up photograph of a black, padded foot fixing device attached to a white plastic foot rest. The device has a silver metal band across the foot and a black strap with a buckle. The brand name "Vey" is visible on the side of the device.	<p><b>Foot fixing</b></p> <ul style="list-style-type: none"><li>- recommended for users that want to secure their feet on their own easily and fast in the foot rests (e. g. for spasticity, palsy,...)</li><li>- fitted with soft and variably positionable cushions for securing the legs without pressure marks</li><li>- cushions can be removed easily for washing and sanitation</li><li>- easy handling of the foot fixing by soft pressure on the cushion (= closing) respectively on the release lever (= opening)</li><li>- only in combination with one of the foot rests (A001-413 / A001-409 / A001-411)</li></ul>	<p><b>A001-777</b></p>
 A photograph showing a person's leg in a grey pant leg being secured by a black leg support device. The device has a padded section and a joint rod. The brand name "Vey" is visible on the side of the device.	<p><b>Leg support</b></p> <ul style="list-style-type: none"><li>- leg support for all users with paralytic symptoms in the legs</li><li>- fitted with soft and variably positionable cushions for securing the legs without pressure marks</li><li>- cushions can be removed easily for washing and sanitation</li><li>- flexible mounted leg support enables natural movement of the ankles</li><li>- cushion and joint rod of the leg support can be removed from the foot rests without tools</li><li>- only in combination with one of the foot rests (A001-413 / A001-409 / A001-411)</li></ul>	<p><b>A001-428</b></p>

## THERA-Trainer Tigo – Handle accessories

	<p><b>Grip adapters</b></p> <ul style="list-style-type: none"> <li>- if the upper body exerciser is used as a safety handle during the leg exercise</li> <li>- easily washable and sanitizable</li> <li>- toolfree assembling to the upper body exerciser</li> <li>- only in combination with long cranks (A001-417)</li> </ul>	<p><b>A001-419</b></p>
	<p><b>Therapy grips</b></p> <ul style="list-style-type: none"> <li>- with three different ergonomic grip positions</li> <li>- easily washable and sanitizable</li> <li>- toolfree assembling to the upper body exerciser</li> <li>- only in combination with one of the cranks (A001-417 / A001-418)</li> </ul>	<p><b>A001-420</b></p>
	<p><b>Arm rests for therapy grips</b></p> <ul style="list-style-type: none"> <li>- for the ergonomic guidance of the forearms</li> <li>- arm rests with soft foamed moldings and elastic velcro straps - very easily washable and sanitizable</li> <li>- can be clipped onto the therapy grips toolfree, quick and easy</li> <li>- only in combination with therapy grips (A001-420)</li> </ul>	<p><b>A003-724</b></p>
	<p><b>Arm rests for tetraplegia</b></p> <ul style="list-style-type: none"> <li>- especially for tetraplegia or "functional hand"</li> <li>- only in combination with one of the cranks (A001-417 / A01-418)</li> </ul>	<p><b>A004-241</b></p>
	<p><b>Arm rests with bar-shaped grips</b></p> <ul style="list-style-type: none"> <li>- for the ergonomic guidance of the forearms</li> <li>- securing of the forearms with velcro straps</li> <li>- the jointed mounting allows a lateral evasive movement of the arms</li> <li>- easily washable and sanitizable</li> <li>- toolfree assembling to the upper body exerciser</li> <li>- only in combination with one of the cranks (A001-417 / A001-418)</li> </ul>	<p><b>A002-757</b></p>

	<p><b>Wristbands for arm rest with bar-shaped grips</b></p> <ul style="list-style-type: none"><li>- for the ergonomic guidance of the forearms</li><li>- recommended for people with strong restrictions of movement or spasticity</li><li>- only in combination with arm rests with bar-shaped grips (A002-757)</li></ul>	<b>A001-427</b>
	<p><b>Arm rest with bar-shaped grips incl. wristbands</b></p> <ul style="list-style-type: none"><li>- for the ergonomic guidance of the forearms</li><li>- recommended for people with strong restrictions of movement or spasticity</li><li>- only in combination with one of the cranks (A001-417 oder A001-418)</li></ul>	<b>A001-426</b>
	<p><b>Wristbands for therapy grips / safety handle</b></p> <ul style="list-style-type: none"><li>- Wristbands for therapy grips / safety handle</li><li>- for the easy securing of paralysed / weak hands on the therapy grips or safety handle/grip adapter</li><li>- only in combination with one of the safety handles (A001-415 / A002-273) or the grip adapters (A001-419) or the therapy grips (A001-420)</li></ul>	<b>A002-440</b>