



DEUTSCHER HEALTHCARE



Pictured is WM1014LXG (G=Graphite)

LEXIA HOME RANGE

DEUTSCHER LEXIA RANGE

A Modern and Stylish range of beds designed to suit your home.

These beds have soft padded fabric covered surrounds and head board to give the bed a homely feel, while having electrically adjustable profiling functions required for healthcare use.

The Lexia Bed range is manufactured in Australia from the highest quality materials to ensure your bed is safe and reliable. All Lexia beds are designed to Australian standard AS/NZS 3200.2.38 and are registered with the TGA.

MODELS-

WM1004LX (G, W, L)- Long single 203 x 92 cm mattress base size

WM1014LX (G, W, L)- King single 203 x 107 cm Mattress base size

STANDARD FABRIC COLOUR OPTIONS (please specify at time of order)-

G: Graphite, W: Wheat, L: Linen





LEXIA BED RANGE

SPECIFICATIONS

Safe working load	190kg (Total weight of patient & accessories)
Height range (Mattress Base)	240 – 740mm
Standard Fabric colours ★	GRAPHITE - LINEN - WHEAT
Wheels	100mm dual castor (brake & directional lock)
Electrics	Dewert - German Engineered
Handset	Coloured key pad, easy to operate and clean
Standard Back Rest	Style 1
Warranty Steel Frame	10 years
Warranty Electrics	3 years
Warranty Bed surround & head board	12 months
Hoist clearance	140mm

OPTIONS

Battery pack up
Underbed Light
Bluetooth controllable handset
Alternative Fabric colours (see swatches)
Style 2 - Bedhead
Style 3 - Bedhead

ACCESSORIES

Bed Sticks x 3 styles
Over bed Tables
Bed side lockers
Mattresses



ABOUT DEUTSCHER HEALTHCARE

Deutscher Healthcare is a manufacturing company based in Ballarat, Victoria and has been in continuous operation for over 110 years. In earlier years, the parent company, HF. Deutscher Pty. Ltd., specialised in farm equipment and later moved into producing a range of mowers and outdoor power equipment, some of which is still manufactured today. The company is still owned and operated by the original family.

All products are designed and manufactured to a high quality, using first grade materials to ensure the product you buy meets Australian standard requirements; is safe, reliable and built to give trouble free service for many years to come.

September 2018