



Peak Care Ellipse Superlite Carbon Fibre Rollator Black

Product info:

SIZE		Small	Medium	Large
	Height - min	76cm	83.5cm	92cm
	Height - max	96cm	96cm	105cm
	Length front wheels to back wheels	66cm	76cm	76cm
	Width - open	60cm	60cm	60cm
	Width - folded	26cm	26cm	26cm
	Seat width	45cm	45.5cm	45.5cm
	Seat depth	22cm	22cm	22cm
	Seat height from ground	45cm	52cm	60cm
	Safe working load	120kgs	120kgs	120kgs
	Front wheel size	20cm (8")	24.2cm (9.5")	24.2cm (9.5")
	Back wheel size	20cm (8")	20cm (8")	20cm (8")
	Product Weight	6.26kg	6.63kg	6.69kg



Materials:

Carbon Fibre

Features:

- AU STD. tested 492484
- Lightweight carbon fibre frame
- Comfortable grip handles
- Tool free adjustable height
- Tool free adjustable brake
- Cane holder
- Strong safe working load

Accessories:

Cup Holder
9346376000574



Oxygen Tank Holder
9346376000581



Product Code:

9346376000635 - Small
9346376000611 - Medium
9346376000628 - Large