OVER BED & CHAIR TABLES

Classic

Our versatile overbed table for everyday use. Multiple base options with adjustable positioning, designed to suit Petite to King Single beds and compact spaces.

Height range: 790-1140mm
Tilt top surface area: 570mm x 390mm
Tilt range: 0-70 degrees
Fixed top surface area: 355mm x 245mm
Max load: 30kg

Colour Options:

Charcoal | White Ash





Standard

A 90mm profile base that moves easily across surfaces, ideal for use on thicker pile carpets.

Dimensions: 668mm L x 391mm W
Height adjust: 790mm - 1140mm
Material: Powder Coated Steel

	Product Code
Charcoal	ACOBTH
White Ash	ACOBTHW



Low Profile

A low-profile base designed to fit easily under ultra low beds with minimal clearance.

Dimensions: 696mm L x 392mm W
Height adjust: 750mm - 1100mm
Material: Powder Coated Steel

	Product Code
Charcoal	ACOBTL
White Ash	ACOBTLW



C - Shape

A C-shaped base that fits around recliners or under beds, allowing closer positioning during use.

Fits Chairs up to 800mm wide
Dimensions: 920mm L x 447mm W
Height adjust: 790mm - 1140mm
Material: Powder Coated Steel

	Product Code
Charcoal	ACOBTC
White Ash	ACOBTCW

OVER BED & CHAIR TABLES

Extra Long

Our extended overbed table for wider beds. Longer reach design improves access across Double, Queen, and Split King surfaces from one side.



Table top: 772mm x 427mm
Fixed table piece: 355mm x 245mm
Height range: 790mm - 1140mm
Tilt range: 0-70 degrees
Internal base clearance: 800mm
Overall base dimensions 920mm x 447mm
Max load: 30kg

Colour Options:

Charcoal	White Ash
ACOBTEC	ACOBTEW



Base Option

C Shape Base
A C-shaped base that fits around recliners or under beds, allowing closer positioning during use.
Internal Clearance (800mm) and External Dimensions: Length 920mm x Width 447mm, height 90mm

Material: Powder Coated Steel

Key Features



- **Extended 772mm Tilting Work Surface**
  – Provides increased reach toward the centre of wider beds for improved access.
- **Designed for Wider Bed Use**
  – Suits Double, Queen, and Split King homecare and adjustable bed setups.
- **C-Shape Base Design**
  – Allows closer positioning around wider bed bases with improved clearance.
- **Split Top Functionality**
  – Tilting main surface with fixed side section for stable, practical use.
- **Height Adjustable Range**
  – Adjusts to suit different bed heights and user positioning needs.
- **Easy-Clean Surface Finish**
  – Smooth surface supports quick cleaning in everyday care environments.

OVER BED & CHAIR TABLES

Over Chair

Our low-height table for seated use. Purpose-built for recliners and day chairs where standard overbed tables sit too high and restrict access.

Table top: 570mm x 390mm
Fixed table piece: 355mm x 245mm
Height range: 550mm - 740mm
Tilt range: 0-70 degrees
Internal base clearance: 800mm
Overall base dimensions 920mm x 447mm
Max load: 30kg

Colour Options:

Charcoal	White Ash
ACOBTELC	ACOBTELW





Base Option

C Shape Base
The C-shape base (920mm x 447mm) with 800mm internal clearance is designed to fit around most lift recliner and day chair bases, enabling closer positioning and improved user comfort.
Unlike traditional overbed tables designed primarily for beds, this model is engineered specifically for seated use — supporting:

Lift recliners | Electric lift chairs | Day chairs | High-back lounge chairs | Seated rehabilitation environments

Key Features



- **Extra Low Height Range (550–740mm)**
  – Designed for seated use where standard tables sit too high and restrict access.
- **Purpose-Built for Chair Use**
  – Engineered for recliners, lift chairs, and seated care environments.
- **C-Shape Base with Internal Clearance**
  – Fits around most recliner and day chair bases for closer positioning.
- **Split Top Functionality**
  – Tilting main surface with fixed side section for stable, practical use.
- **Improved Seated Positioning**
  – Reduces need for forward reach or arm elevation during use.
- **Easy-Clean Surface Finish**
  – Smooth surface supports quick cleaning in everyday care environments.

OVER BED & CHAIR TABLES

Tray

Enhance everyday comfort and independence with our non-slip silicone tray, engineered to sit securely on fixed surfaces, providing added stability while containing spills and keeping small items safely in place.

Fixed table piece: 355mm x 245mm x 25mm



Colour Option:

Charcoal

ACOT



Use

This tray is ideal for keeping drinks steady, meals contained, and small personal items—such as glasses, hearing aids, or utensils—securely in place. The integrated compartments allow for intuitive organisation, reducing clutter while improving accessibility. Perfect for both clinical and home care environments, it supports users in maintaining dignity and independence during daily activities such as eating, reading, or relaxing.

Material: High-quality medical-grade silicone

Key Features



- **Durable Silicone Construction**
  – Flexible, washable material designed for repeated daily use.
- **Four Non-Slip Compartments**
  – Keeps items organised, stable, and easy to access.
- **Spill Containment Design**
  – Raised edges help contain liquids and reduce mess.
- **Secure, Anti-Roll Surface**
  – Prevents items from moving or rolling during use.

OVER BED & CHAIR TABLES

Overview

Explore our Icare overbed & chair tables, offering a stable surface for meals, work, and recreation—featuring easy positioning and an easy-clean, reversible split top with a tilting workspace and fixed side section, available in multiple sizes and shapes for everyday comfort and independence.

Classic

Features

- Choice of low, standard, or C-shape base
- 45 mm or 90 mm underbed clearance
- Compatible with height-adjustable beds

Applications

- Fits Petite, Long Single, and King Single beds
- Supports everyday in-bed positioning
- General-purpose solution for most users
- Fits tighter spaces with limited clearance



Extra Long

Features

- Extended 772mm tilting top for increased reach
- C-shape base for wider bed positioning
- Compatible with height adjustable beds

Applications

- Fits Double, Queen, and Split King beds
- Reaches further toward bed centre
- Supports one-sided access in shared beds
- Ideal for wider adjustable bed setups



Over Chair

Features

- Lowered height range (550–740mm) for seated use
- Designed for close positioning without obstruction
- Fits around most lift recliners and day chairs

Applications

- Supports seated care and rehabilitation environments
- Enables access without forward leaning
- Easy to move away when standing or exiting chair



Tray

Features

- Durable silicone tray with non-slip surface
- Four divided compartments for organisation
- Designed to sit securely on fixed table section

Applications

- Helps prevent drinks from spilling
- Contains spills for easier cleaning
- Prevents items from rolling off the table
- Keeps small items organised and accessible

